

**FIG. 1**




Image Processing Workflow

FIG. 2

35 31 34 32 36 33

1612

QC Workbench \*\*\*\*\* PRODUCTION MODE \*\*\*\*\*

File Edit View

Selected Exports/Imports: Pass Fail Fail/Reason QC User Comment Needs Mask Flood No Mask Not Flood

	Expt Name	Scan Date	Process	Pass/Fail	Reason	IP Fail Flags	IP Fail Description	Chip Type	Intensity pAll	Mean AvgDiff	Row 5/3 Gap/DH	Row 5/3 B-Actn	Inten- sity	Intensity Offset 2	Offset Count	Spk Actn R-Square	S/C
1	08132MV274A21	2002-01-17 1: Archive	pass			0000000 00000100	saty, Uneven; MG_U74Av	4391	1522	1.000	0.750	130	2399	3752	0.965		
2	08132MV274B21	2002-01-17 1: Archive	pass			0010000 00000100	Gid, Intensity; MG_U74Bv	410	531	1.030	0.751	123	3718	4157	0.983		
3	08132MV274C21	2002-01-17 1: Archive	fail	Dim Locally		1011000 01001000	and, Uneven; MG_U74Cv	188	10	0.604	0.774	66	1751	3613	0.978		
4	08153MV274A11	2002-01-17 1: Archive	pass			1000000 00000000	Gid, Uneven; MG_U74Av	313	73	0.731	0.374	201	5105	2825	0.992		
5	08153MV274B11	2002-01-17 1: Archive	pass			0000000 00000000	acts, Uneven; MG_U74Bv	292	35	0.710	0.369	218	8315	3463	0.993		
6	08153MV274C11	2002-01-17 1: Archive	pass			0000000 00000000	MG_U74Cv	222	10	0.717	0.400	203	9887	1527	0.992		
7	08153MV274A21	2002-01-17 1: Archive	pass			0000000 00000100	saty, Uneven; MG_U74Av	553	181	0.728	0.707	153	3547	3471	0.984		
8	08153MV274B21	2002-01-17 1: Archive	pass			0001000 00000000	Gid, Uneven; MG_U74Bv	291	57	0.709	0.661	98	3155	4007	0.981		
9	08153MV274C21	2002-01-17 1: Archive	pass			0000000 00001000	Gid, Intensity; MG_U74Cv	443	49	0.774	0.649	174	3663	3165	0.994		
10	08177MV274A11	2002-01-17 1: Archive	pass			0000000 00100100	All, Intensity; MG_U74Av	431	155	1.029	0.540	273	2498	3118	0.911		
11	08177MV274B11	2002-01-17 1: Archive	pass			0000000 00100000	All; MG_U74Bv	370	79	0.992	0.515	336	3030	4087	0.908		
12	08177MV274C11	2002-01-17 1: Archive	fail	Bright Overall		0010000 00100100	Gid, Intensity; MG_U74Cv	531	40	0.992	0.486	441	3777	2742	0.902		
13	08194MV274A11	2002-01-17 1: Archive	pass			0000000 00100100	saty, Uneven; MG_U74Av	543	61	0.951	0.526	637	4456	2614	0.915		
14	08194MV274B11	2002-01-17 1: Archive	pass			0000000 00000000	Gid, Uneven; MG_U74Bv	382	98	0.704	0.632	204	5591	3474	0.990		
15	08194MV274C11	2002-01-17 1: Archive	pass			0000000 00000000	acts, Uneven; MG_U74Cv	283	37	0.709	0.554	166	4314	4103	0.991		
16	08194MV274A21	2002-01-17 1: Archive	pass			0000000 00001000	saty, Gid; MG_U74Av	339	25	0.693	0.605	223	4537	3206	0.988		
17	08194MV274B21	2002-01-17 1: Archive	pass			0000000 00000000	Acts, Gid; MG_U74Bv	258	21	0.828	0.399	157	5069	2702	0.992		
18	08194MV274C21	2002-01-17 1: Archive	pass			0000000 00000000	Acts, Gid; MG_U74Cv	258	21	0.828	0.399	157	5069	2702	0.992		
19	08197MV274A11	2002-01-17 1: Archive	pass			0000000 00000100	saty, Uneven; MG_U74Av	447	115	0.720	0.602	156	4033	3152	0.989		
20	08197MV274B11	2002-01-17 1: Archive	pass			0000000 00000000	acts, Uneven; MG_U74Bv	369	55	0.796	0.586	146	4325	3743	0.984		
21	08197MV274C11	2002-01-17 1: Archive	fail	Bright Overall		0000000 00000100	saty, Uneven; MG_U74Cv	566	76	0.732	0.521	183	5456	3401	0.988		
22	08198MV274A11	2002-01-17 1: Archive	pass			0000100 00000000	Gid, Uneven; MG_U74Av	367	90	0.590	0.493	132	3416	4071	0.990		
23	08198MV274B11	2002-01-17 1: Archive	pass			0000000 01000000	and, Uneven; MG_U74Bv	318	34	0.706	0.595	336	7193	2979	0.986		
24	08198MV274C11	2002-01-17 1: Archive	pass			0000000 00000100	saty, Uneven; MG_U74Cv	402	36	0.700	0.520	0.420	192	4637	1702	0.992	
25	08863MV274A21	2002-01-17 1: Archive	fail	Haze		0000000 00000000	Gid, Uneven; MG_U74Av	259	10	0.530	0.420	192	4637	1702	0.992		
26	09169MV274A21	2002-01-17 1: Archive	pass	Bright Overall		0000000 00000100	saty, Uneven; MG_U74Av	549	135	0.905	0.548	207	4989	3479	0.988		
27	09169MV274B21	2002-01-17 1: Archive	pass			0001000 00000100	and, Intensity; MG_U74Bv	372	49	0.777	0.546	164	2775	2705	0.987		
28	09169MV274C21	2002-01-17 1: Archive	pass			0000010 01000100	and, Intensity; MG_U74Cv	1766	2	-2.047	-0.483	89	1786	5442	0.916		
29	09250MV274A11	2002-01-17 1: Archive	fail	High Background		0000110 00111100	Saturated, Sca; MG_U74Av	91	0	0.002	-0.754	163	6780	912	0.995		
30	09250MV274B11	2002-01-17 1: Archive	fail	Den Overall		0000010 00000100	saty, Uneven; MG_U74Bv	72	1	2.997	3.475	348	7366	1451	0.988		
31	09250MV274C11	2002-01-17 1: Archive	fail	Den Overall		0000010 00000100	acts, Intensity; MG_U74Cv	120	1	-0.182	1.838	454	9244	1084	0.988		
32	09270MV274A21	2002-01-17 1: Archive	fail	Den Overall		1000000 00000000	Anti-acts, Gid; MG_U74Av	333	103	0.368	0.261	135	2242	3828	0.987		
33	09270MV274B21	2002-01-17 1: Archive	fail	Den Locally		0000000 00000000	MG_U74Bv	333	86	0.402	0.277	185	3974	4027	0.988		
34	09270MV274C21	2002-01-17 1: Archive	pass			0000000 00000100	saty, Uneven; MG_U74Cv	439	52	0.409	0.258	237	3516	3145	0.992		
35	09290MV274A21	2002-01-17 1: Archive	pass			0000000 00000000	Intens; MG_U74Av	557	173	0.699	0.469	131	2729	2886	0.990		

34 of 609

Prod34.Lim All2

Start Microsoft Word QC - Microsoft... Exploring - G... ntelnet - No... John Klein - In... Yahoo! - Micro... QC Workben...

2:23 PM

FIG. 3

Filter

OK Apply Cancel Save Delete

Scan Date:  
Start Date: < 01/14/2002 > End Date: < 01/21/2002 > History

Process Date:  
Start Date: < 01/21/2002 > End Date: < 01/21/2002 > Scanner Setting: ☒ High ☐ Low ☐ Both Needs Mask: ☐ Yes ☒ No

Experiment Name contains:

Lot Number:

Process:

Problem:

Pass/Fail:

Probe Array:

Fail Reason:

Image Processing Parameters

No Metrics No P/F Flags

IP Test	Or	Low Limit	Hi Limit	Pass	Fail	Both	IP Test	Low Limit	Hi Limit	Pass	Fail	Both
<input type="checkbox"/> Any							<input type="checkbox"/> Right Outer Edge					
<input type="checkbox"/> Avg Intensity (All)							<input type="checkbox"/> Top 25% Edge					
<input type="checkbox"/> Outlier Count							<input type="checkbox"/> Bottom 25% Edge					
<input type="checkbox"/> Saturation Count							<input type="checkbox"/> Left 25% Edge					
<input type="checkbox"/> Spikein R-squared							<input type="checkbox"/> Right 25% Edge					
<input type="checkbox"/> Vert 10% peak/med							<input type="checkbox"/> Top 75% Edge					
<input type="checkbox"/> Avg OligoB2 Intens							<input type="checkbox"/> Bottom 75% Edge					
<input type="checkbox"/> Avg Spikein Intens							<input type="checkbox"/> Left 75% Edge					
<input type="checkbox"/> Spikein Intercept							<input type="checkbox"/> Right 75% Edge					
<input type="checkbox"/> Spikein Slope							<input type="checkbox"/> Horiz 25% Max/Min					
<input type="checkbox"/> Negative PM-MM							<input type="checkbox"/> Vert 25% Max/Min					
<input type="checkbox"/> Vert Outlier Var							<input type="checkbox"/> Horiz 75% Max/Min					
<input type="checkbox"/> Horiz Outlier Var							<input type="checkbox"/> Vert 75% Max/Min					
<input type="checkbox"/> Top Outlier Edge							<input type="checkbox"/> Image 5% Intensity					
<input type="checkbox"/> Bottom Outlier Edge							<input type="checkbox"/> 5/3 B-Actin					
<input type="checkbox"/> Left Outlier Edge							<input type="checkbox"/> 5/3 Gapdh					
							<input type="checkbox"/> Mean Avg Diff					

Page 13 Sec 3 13/24 [At 1 Ln 1 Col 1] REC TRK EXT JVR WPH [X] [Y] Yahoo! - Micro QC Workbench ES Filter 227 PM

**FIG. 4**

401				404		405		406
Site	Chip	ChipType	Tissue Type	IP Fail Count	IP Fail Description	Intensity All	Intens/BG	Image 5%
A	1	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	391	5.446	72
A	2	RG U34A	LIVER, NOS	3	All, Artifacts, Uneven	463	6.88	67
A	3	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	529	7.317	72
A	4	RG U34A	LIVER, NOS	1	All	365	5.105	72
A	5	RG U34A	LIVER, NOS	1	All	556	6.864	81
A	6	RG U34A	LIVER, NOS	1	All	469	7.528	62
A	7	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	469	6.312	74
A	8	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	551	8.636	64
A	9	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	369	6.119	60
A	10	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	453	7.906	57
A	11	RG U34A	LIVER, NOS	1	All	337	4.849	70
A	12	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	362	5.3	68
A	13	RG U34A	LIVER, NOS	1	All	333	4.69	71
A	14	RG U34A	LIVER, NOS	5	All, Artifacts, Grid, Uneven	667	9.038	74
A	15	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	321	4.735	68
B	1	RG U34A	LIVER, NOS	0		315	6.389	49
B	2	RG U34A	LIVER, NOS	0		295	4.836	61
B	3	RG U34A	LIVER, NOS	0		315	5.526	57
B	4	RG U34A	LIVER, NOS	0		306	5.741	53
B	5	RG U34A	LIVER, NOS	2	All, Haze Band	294	5.043	58
B	6	RG U34A	LIVER, NOS	0		305	5.865	52
B	7	RG U34A	LIVER, NOS	2	Artifacts, Uneven	419	6.135	68
B	8	RG U34A	LIVER, NOS	1	Artifacts, Uneven	373	6.602	57
B	9	RG U34A	LIVER, NOS	11	Artifacts, Haze Band, Intensity, Uneven	1294	12.686	102
B	10	RG U34A	LIVER, NOS	2	All, Artifacts, Uneven	508	6.773	75
B	11	RG U34A	LIVER, NOS	1	All	294	5.845	50
B	12	RG U34A	LIVER, NOS	2	Artifacts, Uneven	399	7.087	56
B	13	RG U34A	LIVER, NOS	1	Artifacts, Uneven	417	6.915	60
B	14	RG U34A	LIVER, NOS	0		296	6.004	49
B	15	RG U34A	LIVER, NOS	1	All	208	3.014	69

FIG. 4A

Negative PM-MM	Mean AvgDiff	Raw 5/3' GapDH	Raw 5/3' B- ActIn	Intensity SpikesIn	Intensity OligoB2	Li/Wong Outliers	Saturation Count	Vert 10% peak/med	SpikesIn Intercept
46327	118	0.935	0.685	156	11733	2619	0	1.095	3.597
46522	153	0.847	0.721	238	10182	2690	0	1.105	4.146
45037	194	0.942	0.804	324	10683	2666	0	1.075	4.242
47464	107	0.95	0.729	188	11770	2681	0	1.086	4.143
46077	182	0.829	0.741	271	10916	2676	0	1.085	4.23
45626	164	0.897	0.684	427	11569	2856	2	1.073	4.934
45838	170	0.889	0.792	298	7980	2687	0	1.078	3.88
44793	195	0.873	0.684	197	10089	2823	0	1.068	4.244
45615	123	0.965	0.626	286	9582	2615	0	1.054	5.043
44523	175	0.953	0.734	243	9574	2900	0	1.094	3.806
47193	107	1.016	0.687	214	11018	2758	0	1.065	4.309
47243	101	0.909	0.699	188	12237	2620	0	1.084	3.113
46881	105	0.904	0.75	187	10477	2582	0	1.08	3.368
44376	270	0.935	0.75	331	11284	2646	0	1.108	4.146
47765	85	0.944	0.7	178	11579	2574	0	1.127	3.758
46492	110	0.866	0.692	61	2906	2598	0	1.089	4.157
47925	94	0.901	0.711	64	2652	2349	0	1.088	4.214
47385	103	0.852	0.695	70	3113	2436	0	1.166	4.345
46434	102	0.862	0.689	57	2332	2533	0	1.149	4.1
47087	85	0.88	0.723	65	2629	2421	0	1.265	4.269
45981	98	0.861	0.725	73	2835	2511	0	1.152	4.34
47663	130	0.854	0.747	75	3241	2384	0	1.064	4.382
47533	111	0.845	0.746	88	2946	2373	0	1.058	4.623
51528	266	1.012	0.581	104	4521	1761	1	1.357	4.969
46733	174	0.906	0.73	108	4516	2508	0	1.098	4.693
46909	90	0.897	0.688	54	2432	2411	0	1.13	4.056
46670	128	0.813	0.717	72	3478	2483	0	1.071	4.383
47504	118	0.866	0.659	90	4211	2347	0	1.072	4.68
46688	94	0.837	0.678	75	2793	2562	1	1.12	4.44
49044	62	0.681	0.693	48	3336	2404	0	1.046	3.894

FIG. 4B

Spikeln Slope	Vert Outlier Var	Horiz Outlier Var	Top Outlier Edge Ratio	Bottom Outlier Edge Ratio	Left Outlier Edge Ratio	Right Outlier Edge Ratio	25% Top Edge Ratio	25% Bottom Edge Ratio	25% Left Edge Ratio
0.653	22.646	13.817	0.836	0.913	0.801	0.946	1.182	0.947	1
0.64	29.924	17.055	0.764	0.897	1.079	0.928	1.197	0.95	1.022
0.707	24.722	18.077	0.9	0.917	1.139	0.878	1.18	0.948	1.062
0.497	28.449	16.471	0.908	0.824	0.819	0.938	1.137	0.981	1.046
0.699	28.735	23.511	0.81	0.86	1.15	0.884	1.131	0.956	1.043
0.666	30.882	17.744	0.767	0.908	0.952	0.858	1.138	0.947	1.025
0.81	16.837	19.078	0.847	0.963	1.012	0.995	1.167	0.959	1.016
0.716	30.097	23.347	0.774	0.94	0.989	0.957	1.223	0.966	1.047
0.424	29.957	15.748	0.886	0.98	0.883	0.797	1.092	1.029	1.014
0.779	29.609	18.137	0.811	0.811	1.083	0.95	1.207	0.976	1.064
0.501	24.843	19.022	0.84	0.938	0.888	0.92	1.156	0.968	1.014
0.851	24.225	16.097	0.852	0.878	0.759	0.911	1.155	0.969	0.988
0.807	24.343	8.665	0.841	0.78	0.943	1.03	1.145	0.979	1.02
0.823	25.686	23.463	0.83	0.83	1.19	0.743	1.282	0.962	1.101
0.641	23.967	14.159	0.858	0.875	0.902	0.972	1.199	0.927	1.038
0.941	24.777	22.401	0.802	0.802	0.96	0.908	1.14	0.976	0.953
0.93	21.79	17.26	0.837	0.876	1.066	0.883	1.076	1.013	0.991
0.91	23.558	16.515	0.764	0.949	0.994	0.873	1.077	1.023	0.977
0.93	23.146	12.457	0.779	0.832	1.07	0.886	1.068	1.004	1.022
0.889	18.967	17.122	0.838	0.996	1.113	0.919	1.034	0.972	1.028
0.932	25.84	21.033	0.786	0.849	1.052	1.008	1.106	0.998	1.02
0.92	19.927	17.095	0.696	0.912	0.936	0.889	1.058	1.022	0.898
0.863	25.018	18.476	0.823	0.927	1.073	1.016	1.091	1.012	0.947
0.764	19.74	9.226	0.943	1.08	1.144	0.979	0.742	0.568	0.885
0.957	20.442	18.247	0.887	0.923	0.993	0.905	1.153	1.05	1.05
0.915	15.961	12.333	0.813	0.86	1.015	0.923	1.065	0.991	1.018
0.899	26.67	18.798	0.828	0.864	0.981	0.918	1.102	1.017	0.936
0.834	18.036	18.717	0.874	0.96	1.019	0.952	1.114	1.012	0.955
0.887	25.912	17.238	0.746	0.895	1.051	0.955	1.148	0.993	1.037
0.92	15.16	10.043	0.935	0.935	1.141	1.022	1.058	0.988	1.02

FIG. 4C

25% Right Edge Ratio	75% Top Edge Ratio	75% Bottom Edge Ratio	75% Left Edge Ratio	75% Right Edge Ratio	Horiz 25% Max/Min	Vert 25% Max/Min	Horiz 75% Max/Min	Vert 75% Max/Min	Affy Outliers
1.082	1.494	0.888	0.981	1.16	1.197	1.475	1.648	2.554	28
1.108	1.531	0.904	1.019	1.25	1.195	1.535	1.645	2.69	18
1.055	1.541	0.855	1.185	0.967	1.146	1.533	1.365	2.802	10
1.084	1.374	0.944	1.042	1.185	1.146	1.375	1.522	2.287	46
1.057	1.413	0.9	1.112	1.093	1.143	1.445	1.476	2.46	6
1.078	1.441	0.875	1.078	1.105	1.203	1.428	1.526	2.541	45
1.091	1.486	0.93	0.999	1.202	1.179	1.517	1.639	2.651	18
1.102	1.613	0.909	1.088	1.201	1.192	1.629	1.583	3.077	40
1.128	1.309	1.043	1.051	1.071	1.214	1.348	1.439	2.507	63
1.02	1.556	0.895	1.105	1.099	1.139	1.546	1.478	2.946	194
1.108	1.556	0.884	1.031	1.174	1.187	1.359	1.541	2.607	111
1.115	1.476	0.895	0.963	1.136	1.227	1.373	1.668	2.477	118
1.068	1.402	0.929	1.034	1.118	1.122	1.35	1.489	2.347	112
1.018	1.828	0.863	1.276	0.902	1.211	1.69	1.598	3.552	27
1.079	1.612	0.783	1.09	1.139	1.146	1.456	1.453	2.711	28
1.01	1.41	0.941	0.981	1.063	1.176	1.363	1.447	2.318	50
1.003	1.413	0.961	1.012	1.058	1.111	1.308	1.328	2.355	28
1.03	1.437	0.961	1.057	1.034	1.135	1.385	1.364	2.521	46
0.981	1.406	0.95	1.049	1.001	1.109	1.378	1.316	2.372	18
1.002	1.333	0.965	1.169	1.017	1.115	1.445	1.336	2.219	65
1.024	1.442	0.941	1.035	1.125	1.118	1.42	1.447	2.554	38
1.015	1.305	0.947	0.881	1.01	1.237	1.419	1.539	2.317	23
0.928	1.395	0.952	0.922	0.861	1.188	1.41	1.434	2.478	77
1.354	0.999	0.627	0.91	1.321	2.146	2.741	2.074	2.323	575
1.014	1.477	0.979	1.071	1.054	1.161	1.529	1.423	2.745	104
1.005	1.33	0.905	1.087	0.987	1.121	1.38	1.346	2.274	61
0.947	1.42	0.936	0.873	0.926	1.183	1.408	1.529	2.557	33
0.955	1.37	0.928	0.933	0.93	1.189	1.463	1.474	2.403	55
1	1.465	0.924	1.02	1.075	1.116	1.421	1.368	2.58	20
1.032	1.444	1.007	1.106	1.054	1.087	1.148	1.38	2.26	133

FIG. 4D



[illegible]

**FIG. 4E**

00/523499

IP Fail Flags			
00000010	00000000	00000000	00100000
00000010	10000000	00000000	00100000
00000010	00000000	00000000	00100000
00000000	00000000	00000000	00100000
00000000	00000000	00000000	00100000
00000000	00000000	00000000	00100000
00000010	00000000	00000000	00100000
00000010	00000000	00000000	00100000
00000000	00001000	00000000	00100000
00000010	00000000	00000000	00100000
00000000	00000000	00000000	00100000
00000000	00001000	00000000	00100000
00000000	00000000	00000000	00100000
00000010	01010000	01000000	00100000
00000000	00100000	00000000	00100000
00000000	00000000	00000000	00000000
00000000	00000000	00000000	00000000
00000000	00000000	00000000	00000000
00000000	00000000	00000000	00000000
00000000	00000000	00000000	01100000
00000000	00000000	00000000	00000000
00000000	01000100	00000000	00000000
00000000	00000100	00000000	00000000
00000111	10101111	00000000	01000100
00000010	00000000	00000000	00100000
00000000	00000000	00000000	00100000
00000000	01000100	00000000	00000000
00000010	00000000	00000000	00000000
00000000	00000000	00000000	00000000
00000000	00000000	00000000	00100000

FIG. 4F

11/16

Fields	DESCRIPTION
ID	Sequence (Primary Key)
CHIPID	Not used
EXPERIMENTNA	Link to chip table
PROCESSID	From CV PROCESS
PERSON	User or Application
DATETIME	Timestamp
HISTORY	CURRENT or HISTORY
PROBLEMID	From CV PROBLEM (>0 if
FILENAME	Filename from Analysis or

**FIG. 5**

ID	DESCRIPTION
ANALYSIS	Analysis
VALIDATE	Validate
IMAGEPRO	Image
VQC	Visual QC
MASK	Mask
VALIDCHP	Validate
IMPORT	Import
PUBLISH	Publish
ARCHIVE	Archive

**FIG. 6**

ID	DESCRIPTION
0	ok
1	DAT file not found after scan
2	CEL file not found after scan
3	DAT file created without DB entry
4	CHP file is not found
5	CEL file has been modified or
6	Analyzed with incorrect parameters
7	Analyzed without visual QC
8	CEL file created without DB entry
9	CHP file created without DB entry
10	CEL file is older than DAT file
11	CHP file is older than CEL file
12	Failed Visual QC
13	Failed Image Processing

**FIG. 7**

12/16

Fields	DESCRIPTION
CHIPID	IPK
EXPERIMENTNAM	From Affv
PERSON	QCUser
PROBEARRAYTY	From Affv
COMMENTS	
QCDATE	Timestamp
LOTNUMBER	From Affv
PASSFAIL	Set by vac user
DATESTAMP	Current date
FAILREASON	Reason chip failed QC. (no longer same as defect reason)
NEEDSMASK	Flag indicating image needs to be masked (set by vac user from QC)
MASKED	Flag indicating has been masked ('Y') or not (blank) Set from Qualms
IP_FAILFLAGS	25 flag bits. 1=corresponding metric is out of range (failed)
IP_FAILDESCRIPT	Description of defects implied by failed metrics
IP_LIMITSVER	Version number of limits used to compute IP_FAILFLAGS
32 IP Metric columns	IP_INTENSALL, IP_INTENS_SPIKE*, IP_INTENSOLIGOB2*, IP_OUTLIERS, IP_SATUR, IP_SPIKEINR2, IP_VERT10, IP_SPIKEINICPT*, IP_SPIKEINSLOPE*, IP_NEGATIVEPP, IP_VERTOUTVAR, IP_HOROUTVAR, IP_TOPOUTEDGE, IP_BOTTOMOUTEDGE, IP_LEFTOUTEDGE, IP_RIGHTOUTEDGE, IP_TOPEDGE25, IP_BOTTOMEDGE25, IP_LEFTEDGE25, IP_RIGHTEDGE25, IP_TOPEDGE75, IP_BOTTOMEDGE75, IP_LEFTEDGE75, IP_RIGHTEDGE75, IP_HOR25MINMAX, IP_VERT25MINMAX, IP_HOR75MINMAX, IP_VERT75MINMAX, IP_INTENSE5TH, IP_53GAPDH*, IP_53BACTIN*, IP_MEANAVDIFF* *= no limits for these metrics

FIG. 8

Fields	DESCRIPTION
DEFECTID	Sequence (Primary Key)
OLDDEFECTDESCR	For historic reasons -- no longer used
CLASS	Defect type
IMAGE	Not used
CHIPID	FK. Link to ChipDefects PK
QUADRANT	Not used
DEFECTDESCRIPTI	New Description, linked to CV. FAILREASON

**FIG. 9**

Fields	DESCRIPTION
DEFECTID	Sequence (Primary Key)
SHAPE	0=rectangle. 1=ellipse
IMAGE LEFT	Defect location in image coordinates
IMAGE RIGHT	
IMAGE TOP	
IMAGE BOTTOM	
GRID LEFT	Defect location in cel file (grid) coordinates
GRID RIGHT	
GRID TOP	
GRID BOTTOM	

**FIG. 10**


14/16

REASON
Bright Locally
Bright Overall
Cracked
Crop Circle
Dim Locally
Dim Overall
Haze Band
Haze
High Background
Incorrect ProbeArray
Incorrect Scanner Setting
No Sample
Other
Scanner Failure
Snow

**FIG. 11**

Column	Table	DESCRIPTION
LOT RUN ID	AFFX PHYSICAL AR	Lot Number
PROBE ARRAY NA	AFFX PHYSICAL AR	Chip Type – used to update ProbeArrayType in Chip
EXP COMMENT	AFFX ARRAY EXPER	Scanner setting (High/Low)
PROJECT NAME	AFFX SAMPLE	Project name
SCANDATE	CHIP HYB SCAN INF	Scan Date
SCANNER	CHIP HYB SCAN INF	Scanner Name
FLUIDICS	CHIP HYB SCAN INF	Fluidics Name
STATION	CHIP HYB SCAN INF	Fluidics Station
POSITION	CHIP HYB SCAN INF	Fluidics Position

**FIG. 12**

**FIG. 13**




FIG. 14